

Lectins on paraffin

Date: Lectin histochemistry protocol

Materials:

1. Paraffin sections of organs: Mouse spleen is a good control
2. Wash buffer : 0.05 M Tris HCl/ 150 mM NaCl pH 8.0 / 0.1% Tween 20
3. Diluting buffer : wash buffer with 1% BSA: 500 mg /50 ml of wash buffer.
4. Add CaCl₂ (10mM) and MnCl₂ (10 mM) to diluting buffer just before assay

Procedure:

1. De-paraffinize and rehydrate sections on slides: xylene-3 changes, 100% alcohol-3 changes; 95% alcohol – 3 changes, 70% alcohol – 3 changes, TBST washes x3
2. Block endogenous peroxidases: 0.3% H₂O₂ for 30 minutes at room temperature followed by 3 buffer washes, if using HRP labeled secondary
3. Block endogenous biotin:
 - a) Overlay with 0.1% avidin/PBS for 15 minutes; TBST washes x3
 - b) Overlay with 0.01% biotin/PBS for 15 minutes; TBST washes x 3
4. **Primary reagent** incubation:
overlay with negative control or with diluted lectin:
 - a) diluting buffer control slide
 - b) biotinylated SNA 1:500 (Vector Labs B1305, Elderberry bark lectin)
 - c) biotinylated PNA 1:1000 (Vector labs B1075, Peanut Agglutinin)
 - d) biotinylated ECA 1: 2500 (Vector labs B1145, Erythrina Cristagalli)
 - e) Biotinylated GSL1 1:2000 (Vector Labs B1105, Griffonia Simplicifolia lectin)
 - f) biotinylated MAL 2 1:500 (Vector labs B1265, Maackia Amurensis lectin 2)
 - g) biotinylated UEA 1:500 (Vector labs B1065, Ulex Europaeus agglutinin I)
 - h) biotinylated TL 1: 1000 (Vector labs B1175, Lycopersicon Esculentum-tomato lectin)

incubate 30 minutes at room temperature; or one hour or overnight in refrigerator,
wash in 3 changes of wash buffer

6. **Secondary** : Overlay with labeled secondary: incubate for 30 minutes at room temperature;
wash in 3 changes of washing buffer.
 - a. Alkaline phosphatase labeled streptavidin (Jackson labs 016-050-084) 1:100 in diluting buffer
OR with
 - b. HRP-streptavidin (Jackson Immunoresearch) 1:500 in diluting buffer
7. Make fresh **substrate**; overlay with substrate for 3-5 minutes (Vector Blue in 0.1M Tris /levimasole--
Vector labs SK-5300 Alkaline phosphatase substrate III)
or with AEC if using HRP labeled secondary. wash with 3 changes of wash buffer
8. **Nuclear counterstain** in nuclear fast red for 30 minutes (Vector labs ----)
OR Meyer's hematoxylin for 3 minutes; wash with 3 changes of wash buffer
9. **coverslip** with aquamount

Reagents:

1. Bovine Serum Albumin: **Sigma: A4503-50G**
2. **Hydrogen Peroxide to block endogenous peroxidases: Fisher: H325-100**
3. Avidin Biotin Blocking agent Vector labs **SP2001**
4. Alkaline phosphatase labeled streptavidin (**Jackson labs 016-050-084**) 1:100
5. Vector Blue Alkaline phosphatase substrate : Vector labs SK-5300 : made fresh in 0.1M Tris /levimasole--
6. Horse radish Peroxidase labeled streptavidin (**Jackson labs 016-030-084**) 1:500
7. AEC Peroxidase substrate kit Vector labs **SK 4200**

Lectins on paraffin

8. Nuclear counterstain if using AEC substrate for HRP: Mayers' Hematoxylin Sigma Cat No MH532-IL
9. Nuclear counterstain if using Alkaline Phosphatase Vector Blue Substrate: Nuclear Fast Red (Vector labs ----)
10. Aqueous mounting with coverslips:

- SNA**: *Sambucus nigra* agglutinin-- α 2-6 linked sialic acids
- **ECA**: *Erythrina crista-galli* agglutinin-- Gal β 1-4 GlcNAc
- PNA**: Peanut agglutinin-- Gal β 1-3 GalNAc
- MAH**: *Maackia Amurensis* Hemagglutinin-- α 2-3 linked sialic acids
- GSL1**: *Griffonia Simplicifolia* Lectin 1---- α galactose
- TL**: Tomato lectin, *Lycopersicon Esculentum*----poly lactosamines
- L-PHA**: *Phaseolus Vulgaris* Leucoagglutinin----complex structures
- UEA**: *Ulex europaeus* agglutinin---- α linked fucose residues

